

Pseudomonas Gene Deletion Protocol Using pUC18-RedS Lambda Red Recombinase System

Updated on April 6th, 2016

Note: This protocol uses P. aeruginosa PA14 as an example. You need to adapt it based on the Pseudomonas strain you are using.

Day 1

1. Aseptically transfer the filter paper with lyophilized plasmid to a 1.5ml tube. Add 50 μ l of TE (pH 8.0) and incubate at room temperature for at least 20 min but ideally about one hour. You can incubate the paper and TE@37°C with periodic vortexing which may speed the rehydration process. Unless otherwise noted, each filter paper contains approximately 400ng of plasmid.
2. Transform 1 – 2 μ l of the plasmid into your desired high efficiency competent *E. coli* strain and plate on selective media (use 100 μ g/ml ampicillin or carbenicillin depending on the strain) and incubate at the appropriate temperature for your strain.

Day 2

1. Check for transformants. You can verify transformation by PCR (see plasmid map on website) or by plasmid prep and electrophoresis. We normally get hundreds of colonies.

Day 3

1. Prepare a culture of *P. aeruginosa* PA14 carrying pUC18-RedS in 5ml LB + 350 μ g/ml Carb. Shake @37°C overnight (O/N).

Day 4

1. Dilute the O/N culture to OD₆₀₀=0.1 in 5ml LB (1ml of culture is required per transformation).
2. Shake @37°C for 1.5 h (will reach a density (OD₆₀₀) of around 0.4).
3. Add 125 μ l of 20% L-arabinose (final concentration of L-arabinose is 0.5%).
4. Shake @37°C for 1.5 h.
5. Transfer 1 ml of culture into a 1.5 ml tube (Prepare one tube per transformation).
6. Centrifuge at 10,000 x g for 3min at 4°C.
7. Discard supernatant and resuspend bacterial pellet in 1ml of **10% Ice-cold glycerol**.
8. Centrifuge at 10,000 x g for 3min at 4°C.
9. Discard supernatant and resuspend bacterial pellet in 100 μ l **10% Ice-cold glycerol** (These cells are now ready for electroporation).
10. Add at least 1 μ g of DNA fragment consisting of two **~500-bp flanking regions** and antibiotic cassette to the competent cells and incubate 10 min on ice (DNA volume is up to 10 μ l).
11. Transfer competent cells into electroporation cuvette (Bio-Rad #165-2086, 0.2 cm).
12. Electroporate cells with a voltage of 2.5 kV (Ec2 mode on Biorad micropulser).
13. Add 900 μ l of LB into cuvette, resuspend cell well, and transfer contents to a 1.5 ml tube.
14. Shake tube @**30°C** for at least 1h (During this step, prepare a dry LB agar plate(s) containing the appropriate antibiotic).
15. Centrifuge bacterial cells at 10,000 x g for 3min.
16. Discard 800 μ l supernatant.
17. Resuspend pellet by pipetting and spread onto LB agar plate containing appropriate antibiotic.
18. Incubate plate @**30°C** O/N.

Day 3

1. Check for transformant colonies (should appear after 20-24hrs).

Note:

1. We found that a 500bp flanking region gives the best efficiency. Shorter flanking regions can reduce efficiency, but it depends on the region used for recombination.
2. We use 10% ice-cold glycerol instead of 10% sucrose because we found that this has a higher efficiency.
3. After electroporation, we incubate cells at 30°C because the lower temperature may increase the chance of getting colonies.
4. For long-term storage of the plasmid in bacteria, we prefer *E. coli* rather than *P. aeruginosa*, as storage in *P. aeruginosa* can increase the risk in genetic drift of the lambda red system.
5. If you have any additional questions, please contact Dr. Rahme (rahme@research.mgh.harvard.edu).